

DICOM 3.0 Conformance Statement

for

Beacon



Laurel Bridge Software, Inc.
302-453-0222
500 Creek View Rd, Suite 200
Newark, DE 19711
info@laurelbridge.com
www.laurelbridge.com

Document Version: 1.5.0
Document Number: LBDC-000135-010500
Last Saved: 11/30/2022 2:38:00 PM

Copyright 2022, Laurel Bridge Software, Inc. All Rights Reserved

1. Conformance Statement Overview

Beacon is a software application which receives and optionally sends medical instances. It implements the necessary DICOM services to receive those instances from other DICOM entities, and then send those instances to other DICOM entities.

Table 1-1 provides an overview of the DICOM network services supported by Beacon.

Table 1-1
Network Services

SOP Classes	User of Service (SCU)	Provider of Service (SCP)
Transfer		
Computed Radiography Image Storage	Yes	Yes
Digital X-Ray Image Storage - For Presentation	Yes	Yes
Digital X-Ray Image Storage - For Processing	Yes	Yes
Digital Mammography X-Ray Image Storage - For Presentation	Yes	Yes
Digital Mammography X-Ray Image Storage - For Processing	Yes	Yes
Digital Intra-oral X-Ray Image Storage - For Presentation	Yes	Yes
Digital Intra-oral X-Ray Image Storage - For Processing	Yes	Yes
CT Image Storage	Yes	Yes
Enhanced CT Image Storage	Yes	Yes
Ultrasound Multi-frame Image Storage	Yes	Yes
MR Image Storage	Yes	Yes
Enhanced MR Image Storage	Yes	Yes
MR Spectroscopy Storage	Yes	Yes
Nuclear Medicine Image Storage	Yes	Yes
Ultrasound Image Storage	Yes	Yes
Secondary Capture Image Storage	Yes	Yes
Multi-frame Single Bit Secondary Capture Image Storage	Yes	Yes
Multi-frame Grayscale Byte Secondary Capture Image Storage	Yes	Yes
Multi-frame Grayscale Word Secondary Capture Image Storage	Yes	Yes
Multi-frame True Color Secondary Capture Image Storage	Yes	Yes
Standalone Overlay Storage (Retired)	Yes	Yes
Standalone Curve Storage (Retired)	Yes	Yes
Standalone Modality LUT Storage (Retired)	Yes	Yes

SOP Classes	User of Service (SCU)	Provider of Service (SCP)
Standalone VOI LUT Storage (Retired)		
Grayscale Softcopy Presentation State Storage SOP Class	Yes	Yes
Color Softcopy Presentation State Storage SOP Class	Yes	Yes
Pseudo-Color Softcopy Presentation State Storage SOP Class	Yes	Yes
Blending Softcopy Presentation State Storage SOP Class	Yes	Yes
X-Ray Angiographic Image Storage	Yes	Yes
X-Ray Radiofluoroscopic Image Storage	Yes	Yes
Raw Data Storage	Yes	Yes
Spatial Registration Storage	Yes	Yes
Spatial Fiducials Storage	Yes	Yes
Twelve Lead ECG Waveform Storage	Yes	Yes
General ECG Waveform Storage	Yes	Yes
Ambulatory ECG Waveform Storage	Yes	Yes
Hemodynamic Waveform Storage	Yes	Yes
Cardiac Electrophysiology Waveform Storage	Yes	Yes
Basic Voice Audio Waveform Storage	Yes	Yes
VL Endoscopic Image Storage	Yes	Yes
Video Endoscopic Image Storage	Yes	Yes
VL Microscopic Image Storage	Yes	Yes
Video Microscopic Image Storage	Yes	Yes
VL Slide-Coordinates Microscopic Image Storage	Yes	Yes
VL Photographic Image Storage	Yes	Yes
Video Photographic Image Storage	Yes	Yes
Ophthalmic Photography 8 Bit Image Storage	Yes	Yes
Ophthalmic Photography 16 Bit Image Storage	Yes	Yes
Stereometric Relationship Storage	Yes	Yes
Basic Text SR Storage	Yes	Yes
Enhanced SR Storage	Yes	Yes
Comprehensive SR Storage	Yes	Yes
Procedure Log Storage	Yes	Yes
Mammography CAD SR Storage	Yes	Yes
Key Object Selection Document Storage	Yes	Yes
Chest CAD SR Storage	Yes	Yes

SOP Classes	User of Service (SCU)	Provider of Service (SCP)
X-Ray Radiation Dose SR Storage	Yes	Yes
Positron Emission Tomography Image Storage	Yes	Yes
Standalone PET Curve Storage (Retired)	Yes	Yes
RT Image Storage	Yes	Yes
RT Dose Storage	Yes	Yes
RT Structure Set Storage	Yes	Yes
RT Beams Treatment Record Storage	Yes	Yes
RT Plan Storage	Yes	Yes
RT Brachy Treatment Record Storage	Yes	Yes
RT Treatment Summary Record Storage	Yes	Yes
Enhanced XA Image Storage	Yes	Yes
Enhanced XRF Image Storage	Yes	Yes
Encapsulated PDF Storage	Yes	Yes
Enhanced MR Color Image Storage	Yes	Yes
Enhanced US Volume Storage	Yes	Yes
General Audio Waveform Storage	Yes	Yes
Arterial Pulse Waveform Storage	Yes	Yes
Respiratory Waveform Storage	Yes	Yes
XA/XRF Grayscale Softcopy Presentation State Storage	Yes	Yes
Breast Tomosynthesis Image Storage	Yes	Yes
Surface Segmentation Storage	Yes	Yes
Lensometry Measurements Storage	Yes	Yes
Autorefractometry Measurements Storage	Yes	Yes
Keratometry Measurements Storage	Yes	Yes
Subjective Refraction Measurements Storage	Yes	Yes
Visual Acuity Measurements Storage	Yes	Yes
Spectacle Prescription Reports Storage	Yes	Yes
Macular Grid Thickness and Volume Report Storage	Yes	Yes
Colon CAD SR Storage	Yes	Yes
Enhanced PET Image Storage	Yes	Yes
Basic Structured Display Storage	Yes	Yes
Color Palette Storage	Yes	Yes
Intravascular Optical Coherence Tomography Image Storage - For Presentation	Yes	Yes

SOP Classes	User of Service (SCU)	Provider of Service (SCP)
Intravascular Optical Coherence Tomography Image Storage - For Processing	Yes	Yes
VL Whole Slide Microscopy Image Storage	Yes	Yes
Ophthalmic Axial Measurements Storage	Yes	Yes
Intraocular Lens Calculations Storage	Yes	Yes
Ophthalmic Visual Field Static Perimetry Measurements Storage	Yes	Yes
Implantation Plan SR Storage	Yes	Yes
DICOS CT Image Storage	Yes	Yes
DICOS Digital X-Ray Image Storage - For Presentation	Yes	Yes
DICOS Digital X-Ray Image Storage - For Processing	Yes	Yes
DICOS Threat Detection Report Storage	Yes	Yes
Eddy Current Image Storage	Yes	Yes
Eddy Current Multi-frame Image Storage	Yes	Yes
RT Beams Delivery Instruction Storage	Yes	Yes
Generic Implant Template Storage	Yes	Yes
Implant Assembly Template Storage	Yes	Yes
Implant Template Group Storage	Yes	Yes
Workflow Management		
Verification	Yes	Yes

2. Table of Contents

1.	Conformance Statement Overview	2
2.	Table of Contents	6
3.	Introduction	7
3.1.	Revision History	7
3.2.	Audience	7
3.3.	Remarks	7
3.4.	Terms and Definitions	7
3.5.	Basics of DICOM Communication	8
3.6.	Abbreviations	9
3.7.	References	10
4.	Networking	11
4.1.	Implementation Model	11
4.1.1.	Application Data Flow	11
4.1.2.	Functional Definition of AE's	11
4.1.2.1.	Functional Definition of SCU Application Entity (SCU AE)	11
4.1.2.2.	Functional Definition of SCP Application Entity (SCP AE)	11
4.1.3.	Sequencing of Real-World Activities	12
4.1.3.1.	Application Workflow for C-Store	12
4.2.	AE Specifications	12
4.2.1.	SCP Application Entity (SCP AE) Specification	12
4.2.1.1.	SOP Classes	12
4.2.1.2.	Association Policies	17
4.2.1.3.	Association Initiation Policy	18
4.2.1.4.	Association Acceptance Policy	18
4.2.2.	SCU Application Entity (SCU AE) Specification	22
4.2.2.1.	SOP Classes	22
4.2.2.2.	Association Establishment Policies	27
4.2.2.3.	Association Initiation Policy	28
4.2.2.4.	Association Acceptance Policy	30
4.3.	Network Interfaces	30
4.3.1.	Physical Network Interface	30
4.3.2.	Additional Protocols	30
4.3.3.	IPv4 and IPv6 Support	30
4.4.	Configuration	31
4.4.1.	AE Title/Presentation Address Mapping	31
4.4.1.1.	Local AE Titles	31
4.4.2.	Parameters	31
5.	Media Interchange	32
6.	Support of Extended Character Sets	32
7.	Security	33
7.1.	Security Profiles	33
7.1.1.	TLS Security	33
7.2.	Association Level Security	33

3. Introduction

3.1. Revision History

Document Version	Date of Issue	Author	Description
1.0	3 Sep 2019	BAC	Initial creation

3.2. Audience

This document is written for the people that need to understand how Beacon will integrate into their healthcare facility. This includes both those responsible for overall imaging network policy and architecture, as well as integrators who need to have a detailed understanding of the DICOM features of the product. This document contains some basic DICOM definitions so that any reader may understand how this product implements DICOM features. However, integrators are expected to fully understand all the DICOM terminology, how the tables in this document relate to the product's functionality, and how that functionality integrates with other devices that support compatible DICOM features.

3.3. Remarks

The scope of this DICOM Conformance Statement is to facilitate integration between Beacon and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard. DICOM by itself does not guarantee interoperability. The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality. This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of different Conformance Statements is just the first step towards assessing interconnectivity and interoperability between the product and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

3.4. Terms and Definitions

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

Abstract Syntax – the information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples: Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.

Application Entity (AE) – an end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.

Application Entity Title – the externally known name of an Application Entity, used to identify a DICOM application to other DICOM applications on the network.

Application Context – the specification of the type of communication used between Application Entities. Example: DICOM network protocol.

Association – a network communication channel set up between Application Entities.

Attribute – a unit of information in an object definition; a data element identified by a tag. The information may be a complex data structure (Sequence), itself composed of lower level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).

Information Object Definition (IOD) – the specified set of Attributes that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The Attributes may be specified as Mandatory (Type 1), Required but possibly unknown

(Type 2), or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.

Joint Photographic Experts Group (JPEG) – a set of standardized image compression techniques, available for use by DICOM applications.

Media Application Profile – the specification of DICOM information objects and encoding exchanged on removable media (e.g., CDs)

Module – a set of Attributes within an Information Object Definition that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.

Negotiation – first phase of Association establishment that allows Application Entities to agree on the types of data to be exchanged and how that data will be encoded.

Presentation Context – the set of DICOM network services used over an Association, as negotiated between Application Entities; includes Abstract Syntaxes and Transfer Syntaxes.

Protocol Data Unit (PDU) – a packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.

Security Profile – a set of mechanisms, such as encryption, user authentication, or digital signatures, used by an Application Entity to ensure confidentiality, integrity, and/or availability of exchanged DICOM data

Service Class Provider (SCP) – role of an Application Entity that provides a DICOM network service; typically, a server that performs operations requested by another Application Entity (Service Class User). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).

Service Class User (SCU) – role of an Application Entity that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)

Service/Object Pair (SOP) Class – the specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management. **Service/Object Pair (SOP) Instance** – an information object; a specific occurrence of information exchanged in a SOP Class. Examples: a specific x-ray image.

Tag – a 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the “group” and the “element”. If the “group” number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element]

Transfer Syntax – the encoding used for exchange of DICOM information objects and messages. Examples: JPEG compressed (images), little endian explicit value representation.

Unique Identifier (UID) – a globally unique “dotted decimal” string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.

Value Representation (VR) – the format type of an individual DICOM data element, such as text, an integer, a person’s name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.

3.5. Basics of DICOM Communication

This section describes terminology used in this Conformance Statement for the non-specialist. The key terms used in the Conformance Statement are highlighted in *italics* below. This section is not a substitute for training about DICOM, and it makes many simplifications about the meanings of DICOM terms.

Two Application Entities (devices) that want to communicate with each other over a network using DICOM protocol must first agree on several things during an initial network “handshake”. One of the two devices must initiate an Association (a connection to the other device), and ask if specific services, information, and encoding can be supported by the other device (Negotiation).

DICOM specifies a number of network services and types of information objects, each of which is called an Abstract Syntax for the Negotiation. DICOM also specifies a variety of methods for encoding data, denoted Transfer Syntaxes. The Negotiation allows the initiating Application Entity to propose combinations of Abstract Syntax and Transfer Syntax to be used on the Association; these combinations are called Presentation Contexts. The receiving Application Entity accepts the Presentation Contexts it supports.

For each Presentation Context, the Association Negotiation also allows the devices to agree on Roles – which one is the Service Class User (SCU - client) and which is the Service Class Provider (SCP - server). Normally the device initiating the connection is the SCU, i.e., the client system calls the server, but not always.

The Association Negotiation finally enables exchange of maximum network packet (PDU) size, security information, and network service options (called Extended Negotiation information).

The Application Entities, having negotiated the Association parameters, may now commence exchanging data. Common data exchanges include queries for worklists and lists of stored images, transfer of image objects and analyses (structured reports), and sending images to film printers. Each exchangeable unit of data is formatted by the sender in accordance with the appropriate Information Object Definition, and sent using the negotiated Transfer Syntax. There is a Default Transfer Syntax that all systems must accept, but it may not be the most efficient for some use cases. Each transfer is explicitly acknowledged by the receiver with a Response Status indicating success, failure, or that query or retrieve operations are still in process.

Two Application Entities may also communicate with each other by exchanging media (such as a CD-R). Since there is no Association Negotiation possible, they both use a Media Application Profile that specifies “pre-negotiated” exchange media format, Abstract Syntax, and Transfer Syntax.

3.6. Abbreviations

AE	Application Entity
AET	Application Entity Title
CSE	Customer Service Engineer
DHCP	Dynamic Host Configuration Protocol
DICOM	Digital Imaging and Communications in Medicine
DNS	Domain Name System
GSDP	Grayscale Standard Display Function
GSPP	Grayscale Softcopy Presentation State
HIS	Hospital Information System
HL7	Health Level 7 Standard
IHE	Integrating the Healthcare Enterprise
IOD	Information Object Definition
IPv4	Internet Protocol version 4
IPv6	Internet Protocol version 6
ISO	International Organization for Standards
JPEG	Joint Photographic Experts Group
LDAP	Lightweight Directory Access Protocol
LUT	Look-up Table
MPEG	Moving Picture Experts Group
MPPS	Modality Performed Procedure Step
MR	Magnetic Resonance Imaging
MSPS	Modality Scheduled Procedure Step
MTU	Maximum Transmission Unit (IP)
MWL	Modality Worklist
NTP	Network Time Protocol
O	Optional (Key Attribute)
OSI	Open Systems Interconnection
PACS	Picture Archiving and Communication System
PDU	Protocol Data Unit

R	Required (Key Attribute)
RIS	Radiology Information System.
SCP	Service Class Provider
SCU	Service Class User
SOP	Service-Object Pair
SPS	Scheduled Procedure Step
TCP/IP	Transmission Control Protocol/Internet Protocol
U	Unique (Key Attribute)
UL	Upper Layer
VL	Visible Light
VR	Value Representation

3.7. References

NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard, available for free at <https://www.dicomstandard.org>

4. Networking

4.1. Implementation Model

4.1.1. Application Data Flow




Figure 1- Application Data Workflow

4.1.2. Functional Definition of AE's

4.1.2.1. Functional Definition of SCU Application Entity (SCU AE)


The SCU AE can be invoked by the SCP AE's receipt of images and subsequent signing of the study containing the images, triggering the transfer of specific images to a remote destination AE. The SCU AE must be correctly configured with the host and port number of any external DICOM AE's that are to be C-STORE retrieval destinations. The Presentation Contexts to use are determined from the headers of the DICOM files to be transferred, as well as the configuration in Beacon. Some conversion of the DICOM image objects is possible if the original Presentation Context is not supported by the remote destination AE or if a transfer syntax change is preferred.

4.1.2.2. Functional Definition of SCP Application Entity (SCP AE)

The SCP AE waits for another application to connect at the presentation address configured for its Application Entity Title. When another application connects, the SCP AE expects it to be a DICOM application. The SCP AE will accept Associations with Presentation Contexts for SOP Classes of the Transfer, and Workflow Management (Verification only) service classes. Any images received on such Presentation Contexts for Transfer may be added to the Beacon database (depending on configuration).

4.1.3. Sequencing of Real-World Activities

4.1.3.1. Application Workflow for C-Store



4.2. AE Specifications

4.2.1. SCP Application Entity (SCP AE) Specification

4.2.1.1. SOP Classes

Beacon provides Standard Conformance to the following DICOM V3.0 SOP Classes:

Table 4.2.1.1-1
SOP Classes for AE Storage

SOP Classes	SOP Class UID	SCU	SCP
Transfer			
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1	No	Yes
Digital X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.1	No	Yes
Digital X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	No	Yes
Digital Mammography X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.2	No	Yes
Digital Mammography X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	No	Yes
Digital Intra-oral X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.3	No	Yes
Digital Intra-oral X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.3.1	No	Yes
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	No	Yes
Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.1	No	Yes
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	No	Yes
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	No	Yes

SOP Classes	SOP Class UID	SCU	SCP
Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.1	No	Yes
MR Spectroscopy Storage	1.2.840.10008.5.1.4.1.1.4.2	No	Yes
Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.20	No	Yes
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	No	Yes
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	No	Yes
Multi-frame Single Bit Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.1	No	Yes
Multi-frame Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2	No	Yes
Multi-frame Grayscale Word Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.3	No	Yes
Multi-frame True Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.4	No	Yes
Standalone Overlay Storage (Retired)	1.2.840.10008.5.1.4.1.1.8	No	Yes
Standalone Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.9	No	Yes
Standalone Modality LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.10	No	Yes
Standalone VOI LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.11	No	Yes
Grayscale Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.1	No	Yes
Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.2	No	Yes
Pseudo-Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.3	No	Yes
Blending Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.4	No	Yes
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1	No	Yes
X-Ray Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2	No	Yes
Raw Data Storage	1.2.840.10008.5.1.4.1.1.66	No	Yes
Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.1	No	Yes
Spatial Fiducials Storage	1.2.840.10008.5.1.4.1.1.66.2	No	Yes
Twelve Lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	No	Yes
General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2	No	Yes
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3	No	Yes
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1	No	Yes
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1	No	Yes
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1	No	Yes

SOP Classes	SOP Class UID	SCU	SCP
VL Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1	No	Yes
Video Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1.1	No	Yes
VL Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2	No	Yes
Video Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2.1	No	Yes
VL Slide-Coordinates Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3	No	Yes
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4	No	Yes
Video Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4.1	No	Yes
Ophthalmic Photography 8 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1	No	Yes
Ophthalmic Photography 16 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.2	No	Yes
Stereometric Relationship Storage	1.2.840.10008.5.1.4.1.1.77.1.5.3	No	Yes
Basic Text SR Storage	1.2.840.10008.5.1.4.1.1.88.11	No	Yes
Enhanced SR Storage	1.2.840.10008.5.1.4.1.1.88.22	No	Yes
Comprehensive SR Storage	1.2.840.10008.5.1.4.1.1.88.33	No	Yes
Procedure Log Storage	1.2.840.10008.5.1.4.1.1.88.40	No	Yes
Mammography CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.50	No	Yes
Key Object Selection Document Storage	1.2.840.10008.5.1.4.1.1.88.59	No	Yes
Chest CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.65	No	Yes
X-Ray Radiation Dose SR Storage	1.2.840.10008.5.1.4.1.1.88.67	No	Yes
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128	No	Yes
Standalone PET Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.129	No	Yes
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1	No	Yes
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2	No	Yes
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3	No	Yes
RT Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.4	No	Yes
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5	No	Yes
RT Brachy Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.6	No	Yes
RT Treatment Summary Record Storage	1.2.840.10008.5.1.4.1.1.481.7	No	Yes
Enhanced XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1.	No	Yes
Enhanced XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2.1	No	Yes
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	No	Yes
Enhanced MR Color Image Storage	1.2.840.10008.5.1.4.1.1.4.3	No	Yes

SOP Classes	SOP Class UID	SCU	SCP
Enhanced US Volume Storage	1.2.840.10008.5.1.4.1.1.6.2	No	Yes
General Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.2	No	Yes
Arterial Pulse Waveform Storage	1.2.840.10008.5.1.4.1.1.9.5.1	No	Yes
Respiratory Waveform Storage	1.2.840.10008.5.1.4.1.1.9.6.1	No	Yes
XA/XRF Grayscale Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.5	No	Yes
Breast Tomosynthesis Image Storage	1.2.840.10008.5.1.4.1.1.13.1.3	No	Yes
Surface Segmentation Storage	1.2.840.10008.5.1.4.1.1.66.5	No	Yes
Lensometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.1	No	Yes
Autorefractometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.2	No	Yes
Keratometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.3	No	Yes
Subjective Refraction Measurements Storage	1.2.840.10008.5.1.4.1.1.78.4	No	Yes
Visual Acuity Measurements Storage	1.2.840.10008.5.1.4.1.1.78.5	No	Yes
Spectacle Prescription Reports Storage	1.2.840.10008.5.1.4.1.1.78.6	No	Yes
Macular Grid Thickness and Volume Report Storage	1.2.840.10008.5.1.4.1.1.79.1	No	Yes
Colon CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.69	No	Yes
Enhanced PET Image Storage	1.2.840.10008.5.1.4.1.1.130	No	Yes
Basic Structured Display Storage	1.2.840.10008.5.1.4.1.1.131	No	Yes
Color Palette Storage	1.2.840.10008.5.1.4.39.1	No	Yes
Intravascular Optical Coherence Tomography Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.14.1	No	Yes
Intravascular Optical Coherence Tomography Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.14.2	No	Yes
VL Whole Slide Microscopy Image Storage	1.2.840.10008.5.1.4.1.1.77.1.6	No	Yes
Ophthalmic Axial Measurements Storage	1.2.840.10008.5.1.4.1.1.78.7	No	Yes
Intraocular Lens Calculations Storage	1.2.840.10008.5.1.4.1.1.78.8	No	Yes
Ophthalmic Visual Field Static Perimetry Measurements Storage	1.2.840.10008.5.1.4.1.1.80.1	No	Yes
Implantation Plan SR Storage	1.2.840.10008.5.1.4.1.1.88.70	No	Yes
DICOS CT Image Storage	1.2.840.10008.5.1.4.1.1.501.1	No	Yes
DICOS Digital X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.501.2.1	No	Yes

SOP Classes	SOP Class UID	SCU	SCP
DICOS Digital X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.501.2.2	No	Yes
DICOS Threat Detection Report Storage	1.2.840.10008.5.1.4.1.1.501.3	No	Yes
Eddy Current Image Storage	1.2.840.10008.5.1.4.1.1.601.1	No	Yes
Eddy Current Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.601.2	No	Yes
RT Beams Delivery Instruction Storage	1.2.840.10008.5.1.4.34.7	No	Yes
Generic Implant Template Storage	1.2.840.10008.5.1.4.43.1	No	Yes
Implant Assembly Template Storage	1.2.840.10008.5.1.4.44.1	No	Yes
Implant Template Group Storage	1.2.840.10008.5.1.4.45.1	No	Yes
Workflow Management			
Verification	1.2.840.10008.1.1	No	Yes

These are the default SOP Classes supported. By altering the configuration it is possible to support additional or fewer SOP Classes.

4.2.1.1.1. Proposed Presentation Contexts

Beacon is capable of proposing or receiving a presentation context consisting of any SOP class listed in Table 4.2.1.1-1 and any transfer syntax listed in Table 4.2.1.1.1-1 below:

NOTE: Not all SOP classes are inherently compatible with all transfer syntaxes. For example, some of the SOP classes are > 8 bit, and some of the transfer syntaxes are 8 bit only.

Table 4.2.1.1-1

Abstract Syntax Name	Transfer Syntax Name	Transfer Syntax UID	SCU	SCP	Ext. Neg.
*	Explicit VR Little Endian	1.2.840.10008.1.2.1	No	Yes	None
*	Implicit VR Little Endian	1.2.840.10008.1.2	No	Yes	None
*	Explicit VR Big Endian	1.2.840.10008.1.2.2	No	Yes	None
*	RLE Lossless	1.2.840.10008.1.2.5	No	Yes	None
*	JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 [Selection Value 1])	1.2.840.10008.1.2.4.70	No	Yes	None
*	JPEG 2000 (Lossless Only)	1.2.840.10008.1.2.4.90	No	Yes	None
*	JPEG 2000	1.2.840.10008.1.2.4.91	No	Yes	None

Abstract Syntax Name	Transfer Syntax Name	Transfer Syntax UID	SCU	SCP	Ext. Neg.
*	JPEG Lossless, Non-Hierarchical (Process 14)	1.2.840.10008.1.2.4.57	No	Yes	None
*	JPEG Baseline (Process 1)	1.2.840.10008.1.2.4.50	No	Yes	None
*	JPEG Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51	No	Yes	None
*	JPEG-LS Lossless	1.2.840.10008.1.2.4.80	No	Yes	None
*	JPEG-LS Lossy (Near-Lossless)	1.2.840.10008.1.2.4.81	No	Yes	None
*	MPEG2 Main Profile @ Main Level	1.2.840.10008.1.2.4.100	No	Yes	None
*	MPEG2 Main Profile @ High Level	1.2.840.10008.1.2.4.101	No	Yes	None
*	JPEG 2000 Part 2 Multi-component (Lossless Only)	1.2.840.10008.1.2.4.92	No	Yes	None
*	JPEG 2000 Part 2 Multi-component	1.2.840.10008.1.2.4.93	No	Yes	None

* Indicates any SOP class from Table 4.2.1.1-1.

NOTE: A SOP Class can be accepted in the following transfer syntaxes, so long as it is sent out in that same transfer syntax; i.e., Beacon cannot decompress it and recompress it into a different transfer syntax.

Transfer Syntax Name	Transfer Syntax UID
JPEG-LS Lossless	1.2.840.10008.1.2.4.80
JPEG-LS Lossy (Near-Lossless)	1.2.840.10008.1.2.4.81
JPEG 2000 Part 2 Multi-component (Lossless Only)	1.2.840.10008.1.2.4.92
JPEG 2000 Part 2 Multi-component	1.2.840.10008.1.2.4.93
MPEG2 Main Profile @ Main Level	1.2.840.10008.1.2.4.100
MPEG2 Main Profile @ High Level	1.2.840.10008.1.2.4.101

4.2.1.2. Association Policies

4.2.1.2.1. General

The SCP AE will accept Association Requests for the Transfer and Workflow Management (Verification only) services. The DICOM standard Application Context Name for DICOM 3.0 is always accepted and proposed:

**Table 4.2.1.2.1-1
DICOM Application Context For STORAGE-SCP AE**

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.2.1.2.2. Number of Associations

The SCP AE can support multiple simultaneous Associations requested by peer AEs. Each time the SCP AE receives an Association, a child process will be spawned to process the Transfer or Workflow Management (Verification only) service requests. The maximum number of child processes, and thus the maximum number of simultaneous Associations that can be processed, is set by configuration. The default maximum number is 25 in total. This maximum number of simultaneous Associations can be either an absolute number or a maximum number for each requesting external Application Entity. The latter flexibility can be useful if communication with one external AE is unreliable and one does not wish 'hung' connections with this AE to prevent Associations with other client AEs.

**Table 4.2.1.2.1-2
Number Of Simultaneous Associations For SCP AE**

Maximum number of simultaneous Associations requested by peer AEs	25 (Configurable)
---	-------------------

4.2.1.2.3. Asynchronous Nature

The SCP AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2.1.2.1-3
Asynchronous Nature For SCP AE**

Maximum number of outstanding asynchronous transactions	1 (Not configurable)
---	----------------------

4.2.1.2.4. Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2.1.2.1-5
DICOM Implementation Class and Version for SCP AE**

Implementation Class UID	1.2.840.114089.1.0.0.3.3.40
Implementation Version Name	DCF 3.4.56

4.2.1.3. Association Initiation Policy

SCP AE never initiates an Association, therefore also does not propose any Presentation Contexts.


4.2.1.4. Association Acceptance Policy

4.2.1.4.1. Activity – Receive Images and Verification Requests

4.2.1.4.1.1 Description and Sequencing of Activity

The SCP AE accepts Associations only if they have valid Presentation Contexts. If none of the requested Presentation Contexts are accepted then the Association Request itself is rejected. It can be configured to only accept Associations with certain hosts (using TCP/IP address or hostname), ports, and/or Application Entity Titles.

Figure 4.2.1.4.1.1-1



The following sequencing constraints illustrated in Figure 4.2.1.4.1.1-1 apply to the SCP AE for handling transfer requests:

1. Peer AE opens an Association with the SCP AE.
2. Peer AE sends zero or more Composite SOP Instances.
3. Peer AE closes the Association.

The SCP AE may reject Association attempts as shown in the Table below. The Result, Source and Reason/Diag columns represent the values returned in the corresponding fields of an ASSOCIATE-RJ PDU (see PS 3.8, Section 9.3.4). The following abbreviations are used in the Source column:

- a. 1 – DICOM UL service-user
- b. 2 – DICOM UL service-provider (ASCE related function)
- c. 3 – DICOM UL service-provider (Presentation related function)

Table 4.2.1.4.1.1-1

Result	Source	Reason/Diag	Explanation
2 – rejected-transient	c	2 – local-limit-exceeded	The (configurable) maximum number of simultaneous Associations has been reached. An Association request with the same parameters may succeed at a later time.
2 – rejected-transient	c	1 – temporary congestion	No Associations can be accepted at this time due to the real-time requirements of higher priority activities (e.g. during image acquisition no Associations will be accepted) or because insufficient resources are available (e.g. memory, processes, threads). An Association request with the same parameters may succeed at a later time.

Result	Source	Reason/Diag	Explanation
1 – rejected-permanent	a	2 – application-context-name-not-supported	The Association request contained an unsupported Application Context Name. An association request with the same parameters will not succeed at a later time.
1 – rejected-permanent	a	7 – called-AE-title-not-recognized	The Association request contained an unrecognized Called AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association initiator is incorrectly configured and attempts to address the Association acceptor using the wrong AE Title.
1 – rejected-permanent	a	3 – calling-AE-title-not-recognized	The Association request contained an unrecognized Calling AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association acceptor has not been configured to recognize the AE Title of the Association initiator.
1 – rejected-permanent	b	1 – no-reason-given	The Association request could not be parsed. An Association request with the same format will not succeed at a later time.

- **4.2.1.4.1.2 Accepted Presentation Contexts**

The default Behavior of the SCP AE supports all of the presentation contexts listed in Table 4.2.1.1.1-1. The SCP AE can be configured to accept a subset or a superset of these presentation contexts by adding or removing SOP Classes and/or transfer syntaxes.

If multiple Transfer Syntaxes are proposed per Presentation Context then only the most preferable Transfer Syntax is accepted. The order of Transfer Syntax preference for the SCP AE is configurable.

- **4.2.1.4.1.3 SOP Specific Conformance for Verification SOP Class**

The SCP AE provides standard conformance to the Verification SOP Class as an SCP.

- **4.2.1.4.1.4 SOP Specific Conformance for Storage SOP Class**

The associated Activity with the Storage service is the storage of medical image data received over the network on a designated hard disk. The SCP AE will return a failure status if it is unable to store the images on to the hard disk.

The SCP AE does not have any dependencies on the number of Associations used to send images to it. Images belonging to more than one Study or Series can be sent over a single or multiple Associations. Images belonging to a single Study or Series can also be sent over different Associations. There is no limit on either the number of SOP Instances or the maximum amount of total SOP Instance data that can be transferred over a single Association. The SCP AE is configured to retain the original DICOM data in DICOM Part 10 compliant file format. The SCP AE is Level 2 (Full) conformant as a Storage SCP. In addition, all Private and SOP Class Extended Elements are maintained in the DICOM format files. In addition to saving all Elements in files, a subset of the Elements is stored in Beacon' database to support the creation and organization of jobs for the SCU AE.

The behavior for handling duplicate SOP Instances is not configurable. Beacon will accept (assuming the presentation context is acceptable) any duplicate SOP Instances and route them as is.

Table 4.2.1.4.1.4-1

Service Status	Further Meaning	Error Code	Reason
Success	Success	0000	The Composite SOP Instance was successfully received, verified, and stored in the system database.
Refused	Out of Resources	A700	Indicates that there was not enough disk space to store the image.
Error	Data Set does not match SOP Class	A900	Indicates that the Data Set does not encode a valid instance of the SOP Class specified. This status is returned if the DICOM Object stream can be successfully parsed but does not contain values for one or more mandatory Elements of the SOP Class. The SCP AE does not perform a comprehensive check, as it only checks a subset of required Elements. In addition, if the SOP Class is for a type of image but the SOP Instance does not contain values necessary for its display then this status is returned.

Service Status	Further Meaning	Error Code	Reason
Error	Cannot understand	C000	Indicates that the SCP AE cannot parse the Data Set into Elements.
Warning	Coercion of Data Elements	B000	Indicates that one or more Element values were coerced.

NOTE: If a failure condition does occur when handling an Association then all images previously received successfully over the Association are maintained in the Beacon database. No previously successfully received images are discarded. Even if an image is successfully received but an error occurs transmitting the C-STORE Response then this final image is maintained rather than discarded. If the loss of an Association is detected then the Association is closed.

The Behavior of SCP AE during communication failure is summarized in the following table:

Table 4.2.1.4.1.4-2

Exception	Reason
Timeout expiry for an expected DICOM Message Request (DIMSE level timeout). i.e. The SCP AE is waiting for the next C-STORE Request on an open Association but the timer expires.	The Association is aborted by issuing a DICOM A-ABORT. If some Composite SOP Instances have already been successfully received then they are maintained in the database. They are not automatically discarded because of a later failure.
Timeout expiry for an expected DICOM PDU or TCP/IP packet (Low-level timeout). i.e. The SCP AE is waiting for the next C-STORE Data Set PDU but the timer expires.	The Association is aborted by issuing a DICOM A-ABORT. If a C-STORE Data Set has not been fully received then the data already received is discarded. If some Composite SOP Instances have already been successfully received over the Association then they are maintained in the database. They are not automatically discarded because of a later failure.
Association aborted by the SCU or the network layers indicate communication loss (i.e. low-level TCP/IP socket closure)	If some Composite SOP Instances have already been successfully received then they are maintained in the database. They are not automatically discarded because of a later failure.

4.2.2. SCU Application Entity (SCU AE) Specification

4.2.2.1. SOP Classes

Beacon provides Standard Conformance to the following DICOM V3.0 SOP Classes:

Table 4.2.2.1-1
SOP Classes for AE Storage

SOP Classes	SOP Class UID	SCU	SCP
Transfer			
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1	Yes	No
Digital X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.1	Yes	No
Digital X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	Yes	No
Digital Mammography X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.2	Yes	No
Digital Mammography X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	Yes	No

SOP Classes	SOP Class UID	SCU	SCP
Digital Intra-oral X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.3	Yes	No
Digital Intra-oral X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.3.1	Yes	No
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Yes	No
Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.1	Yes	No
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	Yes	No
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	Yes	No
Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.1	Yes	No
MR Spectroscopy Storage	1.2.840.10008.5.1.4.1.1.4.2	Yes	No
Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.20	Yes	No
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	Yes	No
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	Yes	No
Multi-frame Single Bit Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.1	Yes	No
Multi-frame Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2	Yes	No
Multi-frame Grayscale Word Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.3	Yes	No
Multi-frame True Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.4	Yes	No
Standalone Overlay Storage (Retired)	1.2.840.10008.5.1.4.1.1.8	Yes	No
Standalone Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.9	Yes	No
Standalone Modality LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.10	Yes	No
Standalone VOI LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.11	Yes	No
Grayscale Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.1	Yes	No
Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.2	Yes	No
Pseudo-Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.3	Yes	No
Blending Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.4	Yes	No
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1	Yes	No
X-Ray Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2	Yes	No
Raw Data Storage	1.2.840.10008.5.1.4.1.1.66	Yes	No
Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.1	Yes	No

SOP Classes	SOP Class UID	SCU	SCP
Spatial Fiducials Storage	1.2.840.10008.5.1.4.1.1.66.2	Yes	No
Twelve Lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	Yes	No
General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2	Yes	No
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3	Yes	No
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1	Yes	No
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1	Yes	No
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1	Yes	No
VL Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1	Yes	No
Video Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1.1	Yes	No
VL Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2	Yes	No
Video Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2.1	Yes	No
VL Slide-Coordinates Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3	Yes	No
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4	Yes	No
Video Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4.1	Yes	No
Ophthalmic Photography 8 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1	Yes	No
Ophthalmic Photography 16 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.2	Yes	No
Stereometric Relationship Storage	1.2.840.10008.5.1.4.1.1.77.1.5.3	Yes	No
Basic Text SR Storage	1.2.840.10008.5.1.4.1.1.88.11	Yes	No
Enhanced SR Storage	1.2.840.10008.5.1.4.1.1.88.22	Yes	No
Comprehensive SR Storage	1.2.840.10008.5.1.4.1.1.88.33	Yes	No
Procedure Log Storage	1.2.840.10008.5.1.4.1.1.88.40	Yes	No
Mammography CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.50	Yes	No
Key Object Selection Document Storage	1.2.840.10008.5.1.4.1.1.88.59	Yes	No
Chest CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.65	Yes	No
X-Ray Radiation Dose SR Storage	1.2.840.10008.5.1.4.1.1.88.67	Yes	No
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128	Yes	No
Standalone PET Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.129	Yes	No
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1	Yes	No
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2	Yes	No
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3	Yes	No
RT Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.4	Yes	No

SOP Classes	SOP Class UID	SCU	SCP
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5	Yes	No
RT Brachy Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.6	Yes	No
RT Treatment Summary Record Storage	1.2.840.10008.5.1.4.1.1.481.7	Yes	No
Enhanced XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1	Yes	No
Enhanced XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2.1	Yes	No
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Yes	No
Enhanced MR Color Image Storage	1.2.840.10008.5.1.4.1.1.4.3	Yes	No
Enhanced US Volume Storage	1.2.840.10008.5.1.4.1.1.6.2	Yes	No
General Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.2	Yes	No
Arterial Pulse Waveform Storage	1.2.840.10008.5.1.4.1.1.9.5.1	Yes	No
Respiratory Waveform Storage	1.2.840.10008.5.1.4.1.1.9.6.1	Yes	No
XA/XRF Grayscale Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.5	Yes	No
Breast Tomosynthesis Image Storage	1.2.840.10008.5.1.4.1.1.13.1.3	Yes	No
Surface Segmentation Storage	1.2.840.10008.5.1.4.1.1.66.5	Yes	No
Lensometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.1	Yes	No
Autorefractometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.2	Yes	No
Keratometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.3	Yes	No
Subjective Refraction Measurements Storage	1.2.840.10008.5.1.4.1.1.78.4	Yes	No
Visual Acuity Measurements Storage	1.2.840.10008.5.1.4.1.1.78.5	Yes	No
Spectacle Prescription Reports Storage	1.2.840.10008.5.1.4.1.1.78.6	Yes	No
Macular Grid Thickness and Volume Report Storage	1.2.840.10008.5.1.4.1.1.79.1	Yes	No
Colon CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.69	Yes	No
Enhanced PET Image Storage	1.2.840.10008.5.1.4.1.1.130	Yes	No
Basic Structured Display Storage	1.2.840.10008.5.1.4.1.1.131	Yes	No
Color Palette Storage	1.2.840.10008.5.1.4.39.1	Yes	No
Intravascular Optical Coherence Tomography Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.14.1	Yes	No
Intravascular Optical Coherence Tomography Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.14.2	Yes	No
VL Whole Slide Microscopy Image Storage	1.2.840.10008.5.1.4.1.1.77.1.6	Yes	No

SOP Classes	SOP Class UID	SCU	SCP
Ophthalmic Axial Measurements Storage	1.2.840.10008.5.1.4.1.1.78.7	Yes	No
Intraocular Lens Calculations Storage	1.2.840.10008.5.1.4.1.1.78.8	Yes	No
Ophthalmic Visual Field Static Perimetry Measurements Storage	1.2.840.10008.5.1.4.1.1.80.1	Yes	No
Implantation Plan SR Storage	1.2.840.10008.5.1.4.1.1.88.70	Yes	No
DICOS CT Image Storage	1.2.840.10008.5.1.4.1.1.501.1	Yes	No
DICOS Digital X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.501.2.1	Yes	No
DICOS Digital X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.501.2.2	Yes	No
DICOS Threat Detection Report Storage	1.2.840.10008.5.1.4.1.1.501.3	Yes	No
Eddy Current Image Storage	1.2.840.10008.5.1.4.1.1.601.1	Yes	No
Eddy Current Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.601.2	Yes	No
RT Beams Delivery Instruction Storage	1.2.840.10008.5.1.4.34.7	Yes	No
Generic Implant Template Storage	1.2.840.10008.5.1.4.43.1	Yes	No
Implant Assembly Template Storage	1.2.840.10008.5.1.4.44.1	Yes	No
Implant Template Group Storage	1.2.840.10008.5.1.4.45.1	Yes	No
Workflow Management			
Verification	1.2.840.10008.1.1	Yes	No

4.2.2.1.1. Proposed Presentation Contexts

Beacon is capable of proposing or receiving a presentation context consisting of any SOP class listed in Table 4.2.2.1-1 and any transfer syntax listed in Table 4.2.2.1.1-1 below:

NOTE: Not all SOP classes are inherently compatible with all transfer syntaxes. For example, some of the SOP classes are > 8 bit, and some of the transfer syntaxes are 8 bit only.

Table 4.2.2.1.1-1

Abstract Syntax Name	Transfer Syntax Name	Transfer Syntax UID	SCU	SCP
*	Explicit VR Little Endian	1.2.840.10008.1.2.1	Yes	No
*	Implicit VR Little Endian	1.2.840.10008.1.2	Yes	No
*	Explicit VR Big Endian	1.2.840.10008.1.2.2	Yes	No
*	RLE Lossless	1.2.840.10008.1.2.5	Yes	No

Abstract Syntax Name	Transfer Syntax Name	Transfer Syntax UID	SCU	SCP
*	JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 [Selection Value 1])	1.2.840.10008.1.2.4.70	Yes	No
*	JPEG 2000 (Lossless Only)	1.2.840.10008.1.2.4.90	Yes	No
*	JPEG 2000	1.2.840.10008.1.2.4.91	Yes	No
*	JPEG Lossless, Non-Hierarchical (Process 14)	1.2.840.10008.1.2.4.57	Yes	No
*	JPEG Baseline (Process 1)	1.2.840.10008.1.2.4.50	Yes	No
*	JPEG Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51	Yes	No
*	JPEG-LS Lossless	1.2.840.10008.1.2.4.80	Yes	No
*	JPEG-LS Lossy (Near-Lossless)	1.2.840.10008.1.2.4.81	Yes	No
*	MPEG2 Main Profile @ Main Level	1.2.840.10008.1.2.4.100	Yes	No
*	MPEG2 Main Profile @ High Level	1.2.840.10008.1.2.4.101	Yes	No
*	JPEG 2000 Part 2 Multi-component (Lossless Only)	1.2.840.10008.1.2.4.92	Yes	No
*	JPEG 2000 Part 2 Multi-component	1.2.840.10008.1.2.4.93	Yes	No

* Indicates any SOP class from Table 4.2.2.1-1.

4.2.2.2. Association Establishment Policies

4.2.2.2.1. General

The SCU AE can only request the opening of an Association. It cannot accept requests to open Associations from external Application Entities.

The DICOM standard Application Context Name for DICOM is always proposed:

Table 4.2.2.2.1-1
DICOM Application Context For SCU AE

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.2.2.2.2. Number of Associations

The maximum number of simultaneous Associations is configurable. The SCU AE can initiate simultaneous Associations to a given external Destination AE up to the maximum number configured. Each Destination AE can have its own configurable maximum number of simultaneous Associations, but not more than the system-wide maximum number of simultaneous Associations will be honored.

If the first attempt to open an Association fails then the SCU AE will reschedule the task to attempt it again after a configurable time delay. The total number of times to attempt Association establishment is configurable, with the default being 3.

**Table 4.2.2.2-1
Number Of Associations For SCU AE**

Maximum number of simultaneous Associations	10 (Configurable)
---	-------------------

4.2.2.2.3. Asynchronous Nature

The SCU AE does not support asynchronous communication (multiple outstanding transactions over a single Association). All Association requests must be completed and acknowledged before a new operation can be initiated.

**Table 4.2.2.2.3-1
Asynchronous Nature for SCU AE**

Maximum number of outstanding asynchronous transactions	1 (Not Configurable)
---	----------------------

4.2.2.2.4. Implementation Identifying Information

**Table 4.2.2.2.4-1
DICOM Implementation Class and Version For SCU AE**

Implementation Class UID	1.2.840.114089.1.0.0.3.3.40
Implementation Version Name	DCF 3.4.56

Note that the SCU AE and SCP AE use the same Implementation Class UID and Implementation Version Name. This Version Name is updated with each new release of the product software, as the different AE versions are never released independently.


4.2.2.3. Association Initiation Policy

4.2.2.3.1. Activity – Send Images to an External Peer AE

4.2.2.3.1.1 Description and Sequencing of Activity

The SCU AE will initiate a new Association to transmit images. An Association Request is sent to the specified C-STORE Destination AE and upon successful negotiation of the required Presentation Context the image transfer is started. In all cases an attempt will be made to transmit all the indicated images in a single Association, but this may not always be possible. The Association will be released when all the images have been sent. If an error occurs during transmission over an open Association then the image transfer is halted, unless configured to continue.

**Figure 4.2.2.3.1.1-1
Sequencing of Activity – Send Images to an External Peer AE**



The following sequencing constraints illustrated in Figure 4.2.2.3.1.1-1 apply to the SCU AE:

1. SCU AE opens a new Association with the Peer Storage SCP AE.
2. SCU AE sends Composite SOP Instances.
3. SCU AE closes the Association.

4.2.2.3.1.2 SOP Specific Conformance for Verification SOP Class

Standard conformance is provided to the DICOM Verification Service Class as an SCU. The Verification Service as an SCU is actually only supported as a diagnostic service tool for network communication issues.

SOP Specific Conformance for Image SOP Classes

Composite DICOM SOP Instances are maintained as DICOM Part 10 compliant files in the Beacon database. The entire set of tags received with the image will be saved in Beacon; this includes all Private and SOP Extended Elements. When a SOP Instance is selected for export, its content will be exported as it was originally received except for a few possible exceptions. Any of the tags whose values can have been altered due to filtering changes administered on Beacon or changes to the state of the image data due to compression can be altered when the SOP Instance is exported.

The SCU AE will exhibit the following Behavior according to the Status Code value returned in a C-STORE Response from a destination C-STORE SCP:

**Table 4.2.2.3.1.2-1
SCU AE C-STORE Response Status Handling Behavior**

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The SCP has successfully stored the SOP Instance. If all SOP Instances in a store job have status Success then the store job is marked as Sent.
Refused	Out of Resources	A700 – A7FF	This is treated as a failure. The store job will be requeued if it has not reached its max attempts.
Error	Data Set does not match SOP Class	A900 – A9FF	This is treated as a failure. The store job will be requeued if it has not reached its max attempts.
Error	Cannot Understand	C000-CFFF	This is treated as a failure. The store job will be requeued if it has not reached its max attempts.

Service Status	Further Meaning	Error Code	Behavior
Warning	Coercion of Data Elements	B000	Beacon can configurable treat this status as success or failure. If treated as failure the store job will be requeued if it has not reached its max attempts.
Warning	Data Set does not match SOP Class	B007	Beacon can configurable treat this status as success or failure. If treated as failure the store job will be requeued if it has not reached its max attempts.
Warning	Elements Discarded	B006	Beacon can configurable treat this status as success or failure. If treated as failure the store job will be requeued if it has not reached its max attempts.
Warning	Attribute List Error	0107	Beacon can configurable treat this status as success or failure. If treated as failure the store job will be requeued if it has not reached its max attempts.
Warning	Attribute Value Out of Range	0116	Beacon can configurable treat this status as success or failure. If treated as failure the store job will be requeued if it has not reached its max attempts.
*	*	Other	This is treated as a failure. The store job will be requeued if it has not reached its max attempts.

**Table 4.2.2.3.1.2-2
SCU AE Communication Failure Behavior**

Exception	Behavior
Timeout expiry for an expected DICOM Message Response (DIMSE level timeout).	The Association is aborted using a DICOM A-ABORT and the store-job is requeued if it has not reached its max attempts.
Timeout expiry for an expected DICOM PDU or TCP/IP packet (Low-level timeout).	The Association is aborted using a DICOM A-ABORT and the store-job is requeued if it has not reached its max attempts.
Association A-ABORTed by the SCP or the network layers indicate communication loss (i.e. low-level TCP/IP socket closure)	The store-job is requeued if it has not reached its max attempts.

4.2.2.4. Association Acceptance Policy

The SCU AE does not accept Associations.

4.3. Network Interfaces

4.3.1. Physical Network Interface

Beacon supports any network interface that implements the TCP/IP protocol.

4.3.2. Additional Protocols

None.

4.3.3. IPv4 and IPv6 Support

This product only supports IPv4 connections.

4.4. Configuration

4.4.1. AE Title/Presentation Address Mapping

4.4.1.1. Local AE Titles

The specification of AE titles, TCP/IP addresses, and ports is configurable.

Table 4.4.1.1-1

Application Entity	Role	Default AE Title	Default TCP/IP Port
SCU AE	SCU	Beacon	<None>
SCP AE	SCP	CALLED	11112

The SCU and SCP Application Entities can be configured to have the same AE title. The SCP Application Entity has a default Application Entity of CALLED. It can be configured to accept any Called AE Title.

4.4.2. Parameters

Table 4.4.2-1

SCP AE Configuration Parameters

Parameter	Configurable	Default Value
Maximum Send PDU Size	Yes	1048576
Maximum Receive PDU Size	Yes	1048576
Maximum number of simultaneous Associations	Yes	25
Time-out waiting on an open Association for the next Request message (DIMSE Request message, Association Close Request, etc.) (DIMSE timeout)	Yes	1 Hour
Archival time of SOP Instances that have been successfully sent to a Peer AE	Yes	72 Hours

Table 4.4.2-2

SCU AE Configuration Parameters

Parameter	Configurable	Default Value
Maximum Send PDU Size	Yes	1048576
Maximum Receive PDU Size	Yes	1048576
Maximum number of simultaneous Associations	Yes	10
SCU AE time-out waiting for a Response to a DIMSE Request message. (DIMSE timeout)	Yes	1 Hour
Number of times a failed send job to a C-STORE Destination is automatically retried	Yes	2 retries (for a total of 3 attempts)
Amount of time to wait between retries to a C-STORE Destination	Yes	1 minute
Amount of time to wait prior to initial send attempt to a C-STORE Destination	Yes	0 minutes

5. Media Interchange

Beacon does not support Media Interchange.

6. Support of Extended Character Sets

All Beacon DICOM applications support the following:

ISO_IR 100 (ISO 8859-1:1987 Latin Alphabet No. 1 supplementary set) for routing decisions and filtering; other character sets can be stored and forwarded.

7. Security

It is assumed that Beacon is used within a secured environment. It is assumed that a secured environment includes at a minimum:

- a. Firewall or router protections to ensure that only approved external hosts have network access to Beacon.
- b. Firewall or router protections to ensure that Beacon only has network access to approved external hosts and services.
- c. Any communication with external hosts and services outside the locally secured environment use appropriate secure network channels (e.g. such as a Virtual Private Network (VPN)). Alternatively, Beacon can be configured to send and receive DIMSE messages via an encrypted mechanism such as TLS.

Other network security procedures such as automated intrusion detection may be appropriate in some environments. Additional security features may be established by the local security policy and are beyond the scope of this conformance statement.

7.1. Security Profiles

7.1.1. TLS Security

Beacon supports the BCP195 TLS Secure Transport Connection Profile (See DICOM PS3.15 2015c Security and System Management Profiles, Appendix B.9) for authentication and encryption of communication between it and other DICOM clients and servers. Beacon supports TLS version 1.2 as required by this profile, as well as TLS versions 1.0 and 1.1 (if configured by the user).

For secure DICOM communication using the Beacon client, one should select any/all of the TLS options under the Destinations tab; one can also configure additional advanced TLS options in this tab. For secure DICOM communication using the Beacon server, one should select any/all of the Enable TLS options under the System->Listeners->Encrypted Listener tab; one can also configure additional advanced TLS options in this tab.

To properly configure incoming TLS connections, on the System->General tab, the TLS Certificate File Path and certificate Password must first be set. Then on the System->Listeners->Encrypted Listeners tab, the Listener Name and Port must be configured. The Listen Port default for encrypted connections is port 2762, as recommended by the standard (chapter 15 section B.9: The BCP195 TLS Secure Transport Connection Profile); it is the port on which Beacon will receive TLS encrypted communications (e.g., DICOM associations). The TLS Certificate File Path should be set to the file system location of the certificate that Beacon should present for identification to clients. It is suggested that the certificate be a standard PKCS#12 certificate, and it must contain an exportable private key. The Password must be set to the password for the private key in the certificate. Note: Using a certificate format that does not password protect the private key allows the password setting to be ignored.

7.2. Association Level Security

The SCP AE can be configured to check the following values when determining whether to accept Association Open Requests:

- Calling AE Title
- Called AE Title
- Remote IP Address
- Application Context